


Wheelchair Dimensions



No.	Description	Note	Measurement in cm
1	Seat mechanism height	Measure from the floor to the highest point of the armrest.	
2	Length without footrests	If possible, remove the footrests and measure the length to the rearmost point.	
3	Overall length	Measure from the frontmost to the rearmost point of the wheelchair.	
4	Overall height	Measure from the floor to the top edge of the handles.	
5	Folded width	Fold the wheelchair completely and measure from push rim to push rim.	
6	Unfolded width	With the wheelchair unfolded, measure from push rim to push rim.	

Note: If you are unsure, you can also send photos of the wheelchair so that we can check compatibility more quickly.

Wheelchair Information

Brand			
Model			
Removable footrests	☐ yes	☐ no	
Power assist	☐ yes	☐ no	If yes, which one: _____